B.Well

B.Well Limited

758 Great Cambridge Road, Enfield, Middlesex, EN1 3PN
United Kingdom

B.Well Limited

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WA-88

Wrist Watch Blood Pressure Monitor
Instructions Manual

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www.b-well.ru

Wrist Watch Blood Pressure Monitor WA-88

Instruction Manual

Table of contents

1. Introduction	3
1.1. Features	3
1.2. Important information about self-measurement	3
2. Important information on the subject of blood-pressure and its measurement	3
2.1. How does high/low blood-pressure arise?	3
2.2. Which values are normal?	4
2.3. What can be done, if regular high/low values are obtained?	5
3. The various components of the blood-pressure monitor	6
4. Putting the blood-pressure monitor into operation	6
4.1. Inserting the batteries	7
4.2. Reading the set date	7
4.3. Setting the time and date	7
5. Carrying out a measurement	9
5.1. Before the measurement	9
5.2. Common sources of error	9
5.3. Fitting the cuff	10
5.4. Measuring procedure	11
5.5. Discontinuing a measurement	11
5.6. Memory – storage and recall of the measurements	11
6. Error messages/malfunctions	12
7. Care and maintenance, recalibration	13
8. Reference to standards	13
9. Technical specifications	14
10. Guarantee	15

1. Introduction

1.1. Features

The blood-pressure monitor WA-88 (with integrated time/date display) is a fully automatic, digital blood-pressure measuring device for use on the wrist, which enables very fast and reliable measurement of the systolic and diastolic blood-pressure as well as the pulse frequency by way of the oscillometric method of measuring.

The device offers a very high and clinical tested measurement accuracy and has been designed to provide a maximum of user-friendliness.

Before using, please read through this instruction manual carefully and then keep it in a safe place. For further questions on the subject of blood-pressure and its measurement, please contact your doctor.

Attention!

1.2. Important information about self-measurement

- Do not forget: self-measurement means control, not diagnosis or treatment. Unusual values must always be discussed with your doctor. Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.
- The pulse display is not suitable for checking the frequency of heart pacemakers!
- In cases of cardiac irregularity (Arrhythmia), measurements made with this instrument should only be evaluated after consultation with the doctor.

Electromagnetic interference.

The device contains sensitive electronic components (Microcomputer). Therefore, avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave cookers). These can lead to temporary impairment of the measuring accuracy.

2. Important information on the subject of blood-pressure and its measurement

2.1. How does high/low blood-pressure arise?

The level of blood-pressure is determined in a part of the brain, the so-called circulatory centre, and adapted to the respective situation by way of feedback via the nervous system.

To adjust the blood-pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood vessels is altered. The latter is effected by way of fine muscles in the blood-vessel walls.

The level of arterial blood-pressure changes periodically during the heart activity: during the «blood ejection» (Systole) the value is maximal (systolic blood-pressure value), at the end of the heart's «rest period» (Diastole) minimal (diastolic blood-pressure value).

The blood-pressure values must lie within certain normal ranges in order to prevent particular diseases.

2.2. Which values are normal?

Blood pressure is too high if at rest, the diastolic pressure is above 90 mmHg and/or the systolic blood-pressure is over 140 mmHg. In this case, please consult your doctor immediately. Long-term values at this level endanger your health due to the associated advancing damage to the blood vessels in your body.

With blood-pressure values that are too low, i.e. systolic values under 105 mmHg and/or diastolic values under 60 mmHg, likewise, please consult your doctor.

Even with normal blood-pressure values, a regular self-check with your blood-pressure monitor is recommended. In this way you can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, please keep a record of the level of your blood pressure by carrying out regular self-measurements at specific times of the day. Show these values to your doctor. Never use the results of your measurements to alter independently the drug doses prescribed by your doctor.

Table for classifying blood-pressure values (units mmHg) according to World Health Organization:

Range	Systolic Blood-pressure	Diastolic Blood-pressure	Measures
Hypotension	lower than 100	lower than 60	Consult your doctor
Normal range	between 100 and 140	between 60 and 90	Self-check
Mild hypertension	between 140 and 160	between 90 and 100	Consult your doctor
Moderately serious hypertension	between 160 and 180	between 100 and 110	Consult your doctor
Serious hypertension	higher than 180	higher than 110	Consult your doctor immediately
Specific systolic hypertension	higher than 140	lower than 90	Consult your doctor

i Further information

- If your values are mostly standard under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called «labile hypertension». Please consult your doctor if you suspect that this might be the case.
- Correctly measured diastolic blood-pressure values above 120 mmHg require immediate medical treatment.

2.3. What can be done, if regular increased/low values are obtained?

- a) Please consult your doctor.
- b) Increased blood-pressure values (various forms of hypertension) are associated long- and medium term with considerable risks to health. This concerns the arterial blood vessels of your body, which are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can be the result. Furthermore, with long-term continuously increased blood-pressure values, the heart will become structurally damaged.
- c) There are many different causes of the appearance of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organic malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.
- d) There are measures which you can take, not only for reducing a medically established high blood pressure, but also for prevention. These measures are part of your general way of life:

A) Eating habits

- Strive for a normal weight corresponding to your age. Reduce overweight!
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

B) Previous illnesses

Follow consistently any medical instructions for treating previous illness such as:

- Diabetes (Diabetes mellitus)
- Fat metabolism disorder
- Gout

C) Habits


- Give up smoking completely
- Drink only moderate amounts of alcohol
- Restrict your caffeine consumption (Coffee)

D) Physical constitution:

- After a preliminary medical examination, do regular sport.
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your sporting activities. He will advise you regarding the type and extent of types of sport that are possible for you.

3. The various components of the blood-pressure monitor

The illustration shows the blood-pressure monitor, consisting of:


Wrist cuff 13.5–19.5 cm, for wrist circumference of 13.5 to 19.5 cm

4. Putting the blood-pressure monitor into operation

4.1. Inserting the batteries.

After you have unpacked your device, first insert the batteries. The battery compartment is located on the left lateral side of the device (see illustration).

- Remove cover as illustrated.
- Insert the batteries (2 x size AAA 1.5V), thereby observing the indicated polarity.
- If the battery warning appears in the display, the batteries are empty  and must be replaced by new ones.

Attention!

- After the battery warning appears, the device is blocked until the batteries have been replaced.
- Please use «AAA» Long-Life or Alkaline 1.5V Batteries. The use of 1.2V Accumulators is not recommended.
- If the blood-pressure monitor is left unused for long periods, please remove the batteries from the device.

Functional check: Hold on the O/I button  down to test all the display elements. When functioning correctly all segments must appear.

4.2. Reading the set date

Please press the «Time»  key.

4.3. Setting the time and date

This blood-pressure monitor incorporates an integrated clock with date display. This has the advantage, that at each measurement procedure, not only the blood-pressure values are stored, but also the exact moment of the measurement.

After new batteries have been inserted, the clock begins to run from the following setting: 2003-1-01 12:00 o'clock.


You must then re-enter the date and current time. For this, please proceed as follows (Example: Entering 2006-06-15 Time 09:30 o'clock):

- 1) Press the TIME  button for at least 3 seconds.

The display now indicates the set year, during which the last two characters blink.

- 2) The correct year can be entered by pressing the MEMORY  button. (Example: 3 x press)

- 3) Press the TIME  button again. The display now switches to the current date, during which the first character (month) blinks.

4) The corresponding month can now be entered by pressing the MEMORY **M** button. (Example: 5 x presses)

5) Press the TIME **T** button again. The last two characters (day) are now blinking.


6) The corresponding day can now be entered by pressing the MEMORY **M** button. (Example: 29 x presses)

7) Press the TIME **T** button again. The display now switches to the current time, during which the first character (Hour) blinks.

8) The corresponding hour can now be entered by pressing the MEMORY **M** button. (Example: 8 x presses)

9) Press the TIME **T** button again. The last two characters (Minutes) now blink.

10) The exact time can now be entered by pressing the MEMORY **M** button. (Example: 30 x presses)

11) Now after all settings have been made, press the TIME **T** button once again. The date is briefly displayed and then the time. The input is now confirmed and the clock begins to run.

Further Information

With each press of the button   (TIME, MEMORY) one input is made (e.g. switching over from hours to minutes mode, or altering the value by +1).

However, if you keep the respective button depressed, you can switch more quickly to the desired mode, or find the desired value respectively.

5. Carrying out a measurement

5.1. Before the measurement:

- Avoid eating, smoking as well as all forms of exertion directly before the measurement. All these factors influence the measurement result. Try and find time to relax by sitting in an armchair in a quite atmosphere for about ten minutes before the measurement.
- Measure always on the same wrist (normally left).
- Attempt to carry out the measurements regularly at the same time of day, since the blood-pressure changes during the course of the day.

5.2. Common sources of error:

Note:

Comparable blood-pressure measurements always require the same conditions!


These are normally always quiet conditions.

- All efforts by the patient to support the arm can increase the blood-pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the wrist artery lies considerably lower (higher) than the heart, an erroneously higher (lower) blood-pressure will be measured! (Each 15 cm difference in height results in a measurement error of 10 mmHg!)
- A loose cuff causes false measurement values.
- With repeated measurements, blood accumulates in the respective arm, which can lead to false results. Correctly executed blood-pressure measurements should therefore first be repeated after a 5 minute pause or after the arm has been held up in order to allow the accumulated blood to flow away (after at least 3 minutes).

5.3. Fitting the cuff


- a) Remove all eventual objects and jewellery (e.g. wristwatch) from the wrist in question.

 **Draw the cuff over the wrist.**


- b) The distance between the cuff and the hand should be approx. 10 mm.

- c) Secure the cuff with the Velcro fastener, so that it lies comfortably and not too tight, whereby no space should remain between the cuff and the wrist.

- d) Lay the arm on a table, with the palm upwards. Support the arm a little with a rest (cushion), so that the cuff rests at about the same height as the heart. Take care, that the cuff lies free. Remain so for 2 minutes sitting quietly, before beginning with the measurement.

- e) In case the measurement is to be carried out on the right wrist, position the cuff accordingly. However, measurements should always be carried out on the same wrist, so that the values can be compared with each other.

5.4. Measuring procedure

After the cuff has been appropriately positioned and the setting of the pressure pre-selection has been undertaken, the measurement can begin:

- Press O/I  button to put the device into operation.
- Press the START-button, – the pump begins to inflate the cuff. In the display, the increasing cuff-pressure is continually displayed.
- After reaching the inflation pressure, the pump stops and the pressure slowly falls away. The cuff-pressure (large characters) is displayed during the measurement. When the device has detected the pulse, the heart ♥ symbol in the display begins to blink and a peep tone is audible for every pulse beat.
- When the measurement has been concluded, a long peep tone sounds. The measured systolic and diastolic blood-pressure values as well as the pulse frequency are now displayed.

Example (Fig.): Systole 125, Diastole 73, Pulse 75

The measurement results are displayed, until you switch the device off. If no button is pressed for 5 minutes, the device switches automatically off, to save the batteries.

5.5. Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the O/I  button can be pressed at any time. The device then immediately lowers the cuff-pressure automatically.



5.6. Memory – storage and recall of the measurements

The blood-pressure monitor automatically stores each of the last 30 measurement values.

By pressing the MEMORY  button, the last measurement (MR1) as well as the further last 29 measurements (MR2, MR3, ..., MR30) can be displayed one after the other.

(MR1: Values of the last measurement)

(MR2-MR30: Values of the measurement before MR1)

Further information

Measurements should not occur soon after each other, since otherwise, the results will be falsified. Wait therefore for several minutes in a relaxed position, sitting or lying, before you repeat a measurement.

6. Error messages/malfunctions

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed (Example: Error No. 2).

Error No.	Possible cause(s)
«Err 1»	No pulse has been detected.
«Err 2»	Unnatural pressure impulses influence the measurement result. Reason: The arm was moved during the Measurement (Artefact).
«Err 3»	The inflation of the cuff takes too long. The cuff is not correctly seated.

Other possible malfunctions and their elimination

If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

Malfunction	Remedy
The display remains empty when the instrument is switched on although the batteries are in place.	<ol style="list-style-type: none">1. Check batteries for correct polarity and if necessary insert correctly.2. If the display is unusual, re-insert batteries or exchange them.
The device frequently fails to measure the blood pressure values, or the values measured are too low (too high).	<ol style="list-style-type: none">1. Check the positioning of the cuff.2. Measure the blood-pressure again in peace and quiet under observance of the details made under point 5.
Every measurement results in different values, although the device functions normally and normal values are displayed.	<ol style="list-style-type: none">1. Please pay attention to the following information as well as the points listed under «Common sources of error» and repeat the measurement.
The blood pressure values measured differ from those measured by the doctor.	<ol style="list-style-type: none">1. Record the daily development of the values and consult your doctor about them.

7. Care and maintenance, recalibration

- a) Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- b) The cuff contains a sensitive air-tight bubble. Handle this carefully and avoid all types of straining through twisting or buckling.
- c) Clean the device with a soft, dry cloth. Do not use petrol, thinners or similar solvent. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. The cuff must not be washed!
- d) **Do not drop** the instrument or treat it roughly in any way. Avoid strong vibrations.
- e) Never open the device! Otherwise the manufacturer calibration becomes invalid!

Periodical recalibration

Sensitive measuring devices must from time to time be checked for accuracy. We therefore recommend a periodical inspection of the static pressure display **every 1 year**.

Your specialist dealer would be pleased to provide more extensive information about this.

8. Reference to standards

Device standard:

- Device corresponds to the requirements of the European standard for non-invasive blood-pressure monitor
EN1060-1 / 12:95
EN1060-3 / 09:97
DIN 58130, NIBP – clinical investigation
ANSI / AAMI SP10, NIBP – requirements

Electromagn. compatibility:

- Device fulfils the stipulations of the European standard EN 60601-1-2

Clinical testing:

- The clinical performance test was carried out in Germany according to the DIN 58130 / 1997 procedure N6 (sequential).
- The stipulations of the EU-Guidelines 93/42/EWG for Medical Products Class IIa have been fulfilled.

9. Technical specifications

Weight:	145g (with batteries + cuff)
Size:	80 x 70 x 75mm (including cuff)
Storage temperature:	-5 to +50°C
Humidity:	15 to 85% relative humidity maximum
Operation temperature:	10 to 40°C
Display:	LCD-Display (Liquid Crystal Display)
Measuring method:	oscillometric
Measuring range:	
- SYS/DIA	30 to 280 mmHg
- pulse	40 to 200 per minute
Cuff pressure display range:	0 to 299 mmHg
Memory:	storing the last 30 measurements automatically
Measuring resolution:	1 mmHg
Accuracy:	
- pressure	within ± 3 mmHg
- pulse	>5 % of the reading
Power source:	2 dry cells (Batteries) UM-4, size AAA 1.5V
Accessories:	storage case, 2 dry cells (Batteries) UM-4, size AAA, instruction manual

10. Guarantee

The blood-pressure monitor is guaranteed for 3 years from date of purchase. This guarantee includes the instrument. The guarantee does not apply to damage caused by improper handling, accidents, not following the operating instructions or alterations made to the instrument by third parties. The guarantee is only valid upon presentation of the guarantee card filled out by the dealer.

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