---

**B.Well Limited**

758 Great Cambridge Road, Enfield, Middlesex, EN1 3PN  
United Kingdom

**B.Well Limited**

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## WA-77

"Talking" Automatic Blood Pressure Monitor  
Instructions Manual

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CE0044

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# "Talking" Automatic Blood Pressure Monitor WA-77

## Instruction Manual

<b>Table of contents</b>	2
<b>1. Introduction</b>	2
1.1. Features of the WA-77.....	2
1.2. Important information about self-measurement.....	2
<b>2. Important information on the subject of blood-pressure and its measurement</b>	2
2.1. How does high/low blood-pressure arise?.....	2
2.2. Which values are normal?.....	3
2.3. What can be done, if regular high/low values are obtained?.....	4
<b>3. The various components of the blood-pressure monitor</b>	5
<b>4. Putting the blood-pressure monitor into operation</b>	6
4.1. Inserting the batteries.....	6
4.2. Using a mains adapter.....	6
4.3. Tube connection.....	7
<b>5. Carrying out a measurement</b>	7
5.1. Before the measurement.....	7
5.2. Common sources of error.....	7
5.3. Fitting the cuff.....	8
5.4. Measuring procedure.....	9
5.5. Memory – displaying the last measurement.....	9
5.6. Discontinuing a measurement.....	9
<b>6. Error messages/malfunctions</b>	10
<b>7. Care and maintenance, recalibration</b>	11
<b>8. Reference to standards</b>	11
<b>9. Technical specifications</b>	12
<b>10. Guarantee</b>	12

## 1. Introduction

---

### 1.1. Features of the WA-77

Your blood pressure monitor WA-77, is a fully automatic, digital blood pressure measuring device for use on the upper arm. This monitor enables very fast and reliable measurement of systolic and diastolic blood pressure as well your pulse by use of oscillometric method.

This device, equipped with large LCD display & phonetic functions, is very convenient and easy to visualize as well as hear the operational status and cuff pressure during or after a measurement.

This device is easy to use and is proven in clinical studies to provide excellent accuracy.

Please read through this instruction manual carefully before use and then keep it in a safe place. For further questions on the subject of blood pressure and its measurement, please contact your doctor.

### Attention!

### 1.2. Important information about self-measurement

- Do not forget: self-measurement means Control, not diagnosis or treatment. Unusual values must always be discussed with your doctor. Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.
- The pulse display is not suitable for checking the frequency of heart pacemakers!
- In cases of cardiac irregularity (Arrhythmia), measurements made with this instrument should only be evaluated after consultation with the doctor.

### Electromagnetic interference.

The device contains sensitive electronic components (Microcomputer). Therefore, avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave cookers). These can lead to temporary impairment of the measuring accuracy.

## 2. Important information on the subject of blood-pressure and its measurement

### 2.1. How does high/low blood-pressure arise?

The level of blood-pressure is determined in a part of the brain, the so-called circulatory centre, and adapted to the respective situation by way of feedback via the nervous system.

To adjust the blood-pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood vessels is altered. The latter is effected by way of fine muscles in the blood-vessel walls.

The level of arterial blood-pressure changes periodically during the heart activity: during the «blood ejection» (Systole) the value is maximal (systolic blood-pressure value), at the end of the heart's «rest period» (Diastole) minimal (diastolic blood-pressure value).

The blood-pressure values must lie within certain normal ranges in order to prevent particular diseases.

## 2.2. Which values are normal?

Blood pressure is too high if at rest, the diastolic pressure is above 90 mmHg and/or the systolic blood-pressure is over 160 mmHg. In this case, please consult your doctor immediately. Long-term values at this level endanger your health due to the associated advancing damage to the blood vessels in your body.

Should the systolic blood-pressure values lie between 140 mmHg and 160 mmHg and/or the diastolic blood-pressure values lie between 90 mmHg and 95 mmHg, likewise, please consult your doctor. Furthermore, regular self-checks will be necessary.

With blood-pressure values that are too low, i.e. systolic values under 105 mmHg and/or diastolic values under 60 mmHg, likewise, please consult your doctor.

Even with normal blood-pressure values, a regular self-check with your blood-pressure monitor is recommended. In this way you can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, please keep a record of the level of your blood pressure by carrying out regular self-measurements at specific times of the day. Show these values to your doctor. Never use the results of your measurements to alter independently the drug doses prescribed by your doctor.

Table for classifying blood-pressure values (units mmHG) according to World Health Organization:

<b>Range</b>	<b>Systolic Blood-pressure</b>	<b>Diastolic Blood-pressure</b>	<b>Measures</b>
Hypotension	lower than 100	lower than 60	Check with your doctor
Normal range	between 100 and 140	between 60 and 90	Self-check
Mild hypertension	between 140 and 160	between 90 and 100	Consult your doctor
Moderately serious hypertension	between 160 and 180	between 100 and 110	Consult your doctor
Serious hypertension	higher than 180	higher than 110	Consult your doctor immediately

## Further information

- If your values are mostly standard under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called «labile hypertension». Please consult your doctor if you suspect that this might be the case.
- Correctly measured diastolic blood-pressure values above 120 mmHg require immediate medical treatment.

### 2.3. What can be done, if regular increased/low values are obtained?

- a) Please consult your doctor.
- b) Increased blood-pressure values (various forms of hypertension) are associated long- and medium term with considerable risks to health. This concerns the arterial blood vessels of your body, which are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can be the result. Furthermore, with long-term continuously increased blood-pressure values, the heart will become structurally damaged.
- c) There are many different causes of the appearance of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organic malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.
- d) There are measures which you can take, not only for reducing a medically established high blood pressure, but also for prevention. These measures are part of your general way of life:

#### A) Eating habits

- Strive for a normal weight corresponding to your age. Reduce overweight!
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

#### B) Previous illnesses

- Follow consistently any medical instructions for treating previous illness such as:
- Diabetes (Diabetes mellitus)
- Fat metabolism disorder
- Gout

#### C) Habits

- Give up smoking completely
- Drink only moderate amounts of alcohol
- Restrict your caffeine consumption (Coffee)

## D) Physical constitution


- After a preliminary medical examination, do regular sport.
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your sporting activities. He will advise you regarding the type and extent of types of sport that are possible for you.

## 3. The various components of the blood-pressure monitor

---

The illustration shows the blood-pressure monitor WA-77, consisting of:


### a) Measuring unit

### b) Cuff: optional

EA-Cuff A, for arm circumference 22–32 cm

EA-Cuff L, for arm circumference 32–42 cm (available as special accessory)


## 4. Putting the blood-pressure monitor into operation

### 4.1. Inserting the batteries

After you have unpacked your device, first insert the batteries. The battery compartment is located on the underside of the device (see illustration).

- Remove cover as illustrated.
- Insert the batteries (4 x size AA 1,5 V), thereby observing the indicated polarity.
- If the battery warning appears in the display, the batteries are empty and must be replaced by new ones.

### ⚠ Attention!

- After the battery warning appears, the device is blocked until the batteries have been replaced.
- Please use AA Long-Life or Alkaline 1.5 V Batteries. The use of 1.2 V Accumulators is not recommended.
- If the blood-pressure monitor is left unused for long periods, please remove the batteries from the device.


### Functional check

Hold on the O/I button down to test all the display elements. When functioning correctly all segments must appear.

### 4.2. Using a mains adaptor

- It is possible to operate this blood-pressure instrument through a mains adaptor (output 6 V DC / 600 mA with DIN plug). Make certain that you use a mains adaptor AD-1024C manufactured by B.Well or Microlife. The use of any other adaptors can make yours warranty void.

- Plug the DIN plug into the socket at the left side of the instrument.
- Plug the mains adaptor into a 240/110 V power socket.
- Test that power is available by pressing the O/I button.


### ⓘ Note

- No power is taken from the batteries while the mains adaptor is connected to the instrument.
- If the mains voltage is interrupted during the measurement (e.g. by removal of the mains adaptor from the mains socket), the instrument must be reset by removing the DIN plug from its socket.
- Please consult a specialist dealer if you have questions relating to the mains adaptor.

### 4.3. Tube connection

- Insert the cuff tube into the opening provided on the left side of the instrument, as shown in the diagram.

## 5. Carrying out a measurement

---

### 5.1. Before the measurement:

- Avoid eating, smoking as well as all forms of exertion directly before the measurement. All these factors influence the measurement result. Try and find time to relax by sitting in an armchair in a quiet atmosphere for about ten minutes before the measurement.
- Remove any garment that fits closely to your upper arm.
- Measure always on the same arm (normally left).
- Attempt to carry out the measurements regularly at the same time of day, since the blood-pressure changes during the course of the day.

### 5.2. Common sources of error:


#### **i** Note

Comparable blood-pressure measurements always require the same conditions! These are normally always quiet conditions.

- All efforts by the patient to support the arm can increase the blood-pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the arm artery lies considerably lower (higher) than the heart, an erroneously higher (lower) blood-pressure will be measured! (Each 15 cm difference in height results in a measurement error of 10 mmHg!)
- Cuffs that are too narrow or too short result in false measurement values. Selecting the correct cuff is of extraordinary importance. The cuff size is dependent upon the circumference of the arm (measured in the centre). The permissible range is printed on the cuff. If this is not suitable for your use, please contact your dealer. Note: only use clinically approved Original-Cuffs!
- A loose cuff or a sideways protruding air-pocket causes false measurement values.
- With repeated measurements, blood accumulates in the respective arm, which can lead to false results. Correctly executed blood-pressure measurements should therefore first be repeated after a 1 minute pause.

### 5.3. Fitting the cuff

a) Pass the end of the cuff (with sewn-in rubber stopper) through the metal stirrup so that a loop is formed. The Velcro closer must be facing outwards. (Ignore this step if the cuff has already been prepared.)


b) Push the cuff over the left upper arm so that the tube points in the direction of the lower arm.


c) Lay the cuff on the arm as illustrated. Make certain that the lower edge of the cuff lies approximately 2 to 3 cm above the elbow and that the rubber tube leaves the cuff on the inner side of the arm. Important! The mark (ca. 3 cm long bar) must lie exactly over the artery which runs down the inner side of the arm.


d) Tighten the free end of the cuff and close the cuff with the closer.

e) There must be no free space between the arm and the cuff as this would influence the result. Clothing must not restrict the arm. Any piece of clothing which does (e.g. a pullover) must be taken off.

f) Secure the cuff with the Velcro closer in such a way that it lies comfortably and is not too tight. Lay the arm on the table (palm upwards) so that the cuff is at the same height as the heart. Make sure that the tube is not kinked.

g) Remain seated quietly for two minutes before you begin the measurement.

#### Comment

If it is not possible to fit the cuff to the left arm, it can also be placed on the right one. However all measurements should be made using the same arm.

## 5.4. Measuring procedure

After the cuff has been appropriately positioned the measurement can begin:

a) Press O/I + Memory button. The device reminds you that "Attention! Measuring starts. Do not move and talk". After which pump begins to inflate the cuff. Rising pressure in the cuff is now shown in the display.

b) After the suitable inflation pressure has been reached, pump stops and pressure gradually falls. The cuff pressures are displayed. In case that the inflation pressure is not sufficient, the monitor automatically re-inflates to a higher level.

c) When the instrument detects a pulse, the heart symbol in the display starts to flash and a beep is heard for every heartbeat.

d) A longer beep is sounded when the measurement has been completed. Then the systolic, diastolic blood pressures and pulse rate now appear in the display and device automatically voices out your measured values, by 4 times.

e) The measurement readings remain on the display until you switch off the instrument. If no button is pressed for a period of 3 minutes the instrument switches itself off in order to preserve the batteries.

### **Note**

If it is expected to take measurements without voice, adjust volume to the smallest by rotating the controller on the right side accordingly, then no voice (except beeper) will be heard in the whole measuring process.


## 5.5. Memory – displaying the last measurement

The measured results are stored in the instrument until a new measurement is carried out or the batteries are removed. You can call up the stored values at any time by pressing the O/I + Memory button for at least 3 seconds, the device then automatically shows as well reads out by 4 times the memory reading.

## 5.6. Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g the patient feels unwell), the O/I + Memory button can be pressed at any time. The device then immediately and automatically lowers cuff pressure, at this time you hear "hope you health and a good mood", press O/I + Memory button immediately, it stops.

## 6. Error messages/malfunctions

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed, the device automatically reminds you "measuring failed, please try it again". (Example Error no. 1)

Error No.	Possible cause(s)
«Err 1»	The systolic pressure was determined but afterwards the cuff pressure fell below 20 mmHg. This situation can occur, for example, if the tube becomes unplugged after the systolic blood pressure has been measured. Further possible cause: No pulse has been detected.3 sec.
«Err 2»	Unnatural pressure impulses influence the measurement result. Reason: The arm was moved during the measurement (Artefact).
«Err 3»	Inflation of the cuff takes too long. The cuff is not correctly seated or the hose connection is not sealed.
«Err 5»	The measured reading indicated an unacceptable difference between systolic and diastolic pressures Take another reading following directions carefully. Contact your doctor if you continue to get unusual readings.
HI	The cuff pressure is too high. * The measured pulse is higher than 200 beats/min.
LO	The measured pulse is lower than 40 beats/min. Repeat the measurement.*

\*If this or any other problem occurs repeatedly, please consult your doctor.

### Other possible malfunctions and their elimination

If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

Malfunction	Remedy
The display remains empty when the instrument is switched on although the batteries are in place.	1. Check batteries for correct polarity and if necessary insert correctly. 2. If the display is unusual, re-insert batteries or exchange them.
The device frequently fails to measure the blood pressure values, or the values measured are too low (too high).	1. Check the positioning of the cuff. 2. Measure the blood-pressure again in peace and quiet under observance of the details made under point 5.

Every measurement results in different values, although the device functions normally and normal values are displayed.	1. Please pay attention to the following information as well as the points listed under «Common sources of error» and repeat the measurement.
The blood pressure values measured differ from those measured by the doctor.	1. Record the daily development of the values and consult your doctor about them.

## 7. Care and maintenance, recalibration

a) Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.

b) The cuff contains a sensitive air-tight bubble. Handle this carefully and avoid all types of straining through twisting or buckling.

c) Clean the device with a soft, dry cloth. Do not use petrol, thinners or similar solvent. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. The cuff must not be washed!

d) Handle the tube carefully. Do not pull on it. Likewise, this is not to be buckled or laid over sharp edges.

e) Do not drop the instrument or treat it roughly in any way. Avoid strong vibrations.

f) Never open the device! Otherwise the manufacturer calibration becomes invalid!

### Periodical recalibration

Sensitive measuring devices must from time to time be checked for accuracy. We therefore recommend a periodical inspection of the static pressure display every 1 year.

Your specialist dealer would be pleased to provide more extensive information about this.

## 8. Reference to standards

### Device standard:

Device corresponds to the requirements of the European standard for non-invasive blood-pressure monitor

EN1060-1 / 12:95

EN1060-3 / 09:97

DIN 58130, NIBP – clinical investigation

ANSI / AAMI SP10, NIBP – requirements

### **Electromagn. compatibility:**

Device fulfils the stipulations of the European standard EN 60601-1-2

### **Clinical testing:**

The clinical performance test was carried out in Germany according to the DIN 58130 / 1997 procedure N6 (sequential).

The stipulations of the EU-Guidelines 93/42/EWG for Medical Products Class IIa have been fulfilled.

## **9. Technical specifications**

---

Weight:	560 g (with batteries + cuff)
Size:	140*110*65 mm
Storage temperature:	-5 to +50°C
Humidity:	15 to 85 % relative humidity maximum
Operation temperature:	10 to 40°C
Display:	LCD-Display (Liquid Crystal Display)
Measuring method:	oscillometric
Measuring range:	
SYS/DIA	30 to 280 mmHg
Pulse:	40 to 200 per minute
Cuff pressure display range:	0 to 299 mmHg
Memory:	storing the last measurement
Measuring resolution:	1 mmHg
Accuracy:	
- pressure	within $\pm 3$ mmHg
- pulse	$\pm 5$ % of the reading
Power source:	4 dry cells (Batteries), size AA 1.5V
AC Adapter (optional)	AD-1024C
Accessories:	cuff «EA-Cuff A» 22-32 cm, 4 dry cells (Batteries), size AA, instruction manual, warranty card.

## **10. Guarantee**

---

B.Well warrants to the original retail purchaser that this product is free from defects in material and workmanship for a period of 3 years from the date of purchase. This warranty includes the instrument. The warranty does not apply to damage caused by improper handling, accidents, not following the operating instructions or alterations made to the instrument by third parties.

The warranty is only valid upon presentation of the warranty card filled out by the dealer or by the service center.

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