

injeTAK™

Precision Cystoscopic Injection Needle



Always the right option



InjeTAK Precision Cystoscopic Injection Needle




Adjustable multiple sided needle point
Precision puncture ensures minimal tissue damage.



25G needle with 24G performance
Minimizes tissue damage and pain while maintaining flow during procedure.



Black indicator tip
Ensures puncture depth and safe penetration.



Unique needle design
Maintains structural integrity throughout the entire procedure.





Easy to use
Ergonomic design allows for easy adjustment and single handed operation.



0-5mm Setting
Protects the needle tip and cystoscope when not in use. Adjustable in length for precision accuracy.

Part #	Length	Depth
DIS196	50cm	0,2,3,5mm
DIS198	35cm	



USA:
Tel.: +1 802 857 1300
800 522 6743
Fax: +1 802 878 1122

CANADA:
Tel.: +1 450 671 5901
888 522 6743
Fax: +1 450 671 7182

UK:
Tel.: +44 1275 372332
Fax: +44 1275 375145

FRANCE:
Tel.: +33 1 43 41 34 22
Fax: +33 1 43 41 38 69

INTERNATIONAL:
Tel.: +1 905 612 1170
Fax: +1 905 612 9731

BRAZIL:
Tel.: +55 11 4134 8480
Fax: +55 11 4208 3837

BELGIUM:
Tel.: +32 3 780 17 38
Fax: +32 3 780 17 39

CHINA:
Tel.: +86 10 8489 8891
Fax: +86 10 8489 8511



www.laborie.com